Cochlear™ Baha® Remote Control User Manual

Manufacturer:
Cochlear Bone Anchored Solutions AB Konstruktionsvägen 14, SE - 435 33 Mölnlycke, Sweden
Tel: +46 31 792 44 00 Fax: +46 31 792 46 95

Regional offices:
Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University NSW 2109, Australia Tel: +61 2 9428 6355 Fax: +61 2 9428 6352
Cochlear Americas 13059 E Peakview Avenue, Centennial, CO 80111, USA
Tel: +1 303 790 9010 Fax: +1 303 792 9025
Cochlear AG EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland
Tel: +41 61 205 0404 Fax: +41 61 205 0405

www.cochlear.com

94498

614040-00

Cochlear, Hear now. And always, Hybrid, Nucleus, SmartSound, the elliptical logo, Nucleus in Chinese characters and Codacs, are either trademarks or registered trademarks of Cochlear Limited. Ardiun, Baha, Baha Divino, Baha Intenso, Baha PureSound, Baha SoftWear, DermatLock and Vistafix are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB. Cochlear, Hear Now. And always and the elliptical logo are either trademarks or registered trademarks of Cochlear Limited. © Cochlear Bone Anchored Solutions AB 2013. All rights reserved. JUN13.




Hear now. And always

Cochlear®

Welcome

Congratulations on your purchase of a Cochlear™ Baha® Remote Control.

Key to device *See figure 1*

1. Display
2. Key Lock
3. TV button
4. Home button
5. Increase volume
6. Program change
7. Activate right hearing instrument only
8. Decrease volume
9. Activate left hearing instrument only
10. Pairing button
11. AAAA battery

Note: The instructions in this manual describe operation of the Cochlear Baha Remote Control with one Cochlear sound processor. If you are a bilateral recipient (using two Cochlear sound processors), note that one Cochlear Baha Remote Control can be used with two Cochlear sound processors. If this is the case, the same instructions apply to both sound processors.

Introduction

The Cochlear Baha Remote Control will allow you to operate your Cochlear Baha wireless sound processors. It is both an elegant and user-friendly Remote Control as well as a status indicator that gives you access to other wireless functions on your sound processor.

If you are using two wireless Baha sound processors, the Cochlear Baha Remote Control can be used to adjust both of your sound processors simultaneously. Or you can control both sound processors individually, allowing you to perfectly fine tune your sound processor to any specific hearing situation.

Please read this manual carefully to fully benefit from the Cochlear Baha Remote Control. Ask your hearing care professional if you have any questions.

For further information you can also visit the www.cochlear.com web site.

For your safety, carefully read the chapter "Important Information".

Warranty

The warranty does not cover defects or damage arising from, associated with, or related to the use of this product with any non-Cochlear processing unit and/or any non-Cochlear implant. See "Cochlear Global Limited Warranty card" for more details.

Key to symbols

The following symbols will be used throughout these documents.

	"Consult instructions for use" or "Operating instructions"
	"Caution" or "Caution, consult accompanying documents"
	CE-mark
	Serial Number
	Refer to instructions/booklet
	ACMA (Australian Communications and Media Authority) compliant
	Please ask your hearing care professional concerning disposal of this device.

Getting Started

Changing program

See figure 2

The "P" button (no. 6 in figure 1) is used to change the programs that your hearing care professional has programmed into your sound processor. The name of the program displayed can be configured by your hearing care professional. Each time you press the "P" button the program advances one step, returning to program 1 at the end of each cycle (e.g., 1, 2, 3, 1). Each program change is accompanied by beep tones in both sound processors and the program currently in use is shown in the display.

Adjusting volume

Adjust volume on two sound processors

See figure 3

The "+" and "-" keys are used to adjust volume (nos. 5 and 8 in figure 1). Volume may be adjusted on both sound processors simultaneously or on each side independently. By default, both sound processors are "active," which is indicated in the display by two arrowheads (< and >) located next to the volume bars. Both sound processors can simultaneously be made louder by pressing the "+" key or softer by pressing the "-" key. Each adjustment is accompanied by beep tones in the sound processors as well as a visual change in the volume bars shown in the display.

Adjusting volume on the left sound processor only

See figure 4

To adjust volume on the left side only, first press the "<" key on the keypad (no. 9 in figure 1). Doing so makes the left sound processor "active" and the right sound processor "inactive". This is indicated by a left arrowhead (<) in the display, with the right arrowhead turned off. Press "+" or "-" to adjust volume on the left sound processor only.

Adjusting volume on the right sound processor only

See figure 5

To adjust volume on the right side only, first press the ">" key on the keypad (no. 7 in figure 1). Doing so makes the right side "active" and the left side "inactive". This is indicated by a right arrowhead (>) in the display with the left arrowhead turned off. Press "+" or "-" to adjust volume on the right side only.

To return to simultaneous adjustment mode, press either the < or > keys on the keypad (nos. 7 or 9 in figure 1) until both arrowheads appear in the display, indicating that both sound processors are again "active".

Return to Home settings

Home settings, i.e., those that apply whenever the sound processor is first turned on, are activated by pressing the Home button (no. 4 in figure 1). Pressing the Home button will:

1. Activate program 1 in your sound processor
2. Adjust your sound processor to the default volume settings
3. Set the volume to simultaneous adjustment mode

The home button will also terminate any streamed audio signals from Cochlear TV Streamer.

Key lock

See figure 6

A two-position key lock switch, located on the right side of the Remote Control, helps prevent accidental activation. This feature is handy when the device is kept in a pocket or purse.

When in the locked position all keys/buttons are non-functional. When switched to the non-locked position normal functionality is restored.

When the key lock is activated, the key lock symbol will appear on the screen when any button is pushed. *See figure 7*

Note: When key lock is activated all wireless Remote Control operation is disabled. Activate the key lock during airplane flights.

Low battery indicator

See figure **8**

In the event that the Remote Control's batteries are weak, this will be indicated by a low-battery icon situated above the program name. The low-battery icon and the program that your sound processor is in will alternate on the display.

Sound processor low battery

When the sound processor's battery is almost flat, you may still be able to use it, but the Remote Control functionality may be disabled. Should this happen, we recommend you replace the battery in your sound processor.

Operating range

See figure **9**

The operating range between the Remote Control and the sound processor is approximately 1 meter (40"). If the current status is not displayed on the Remote Control, hold Remote Control in front of your chest as shown in the illustration.

If the Remote Control cannot determine the sound processor's status, search icons will appear on its display. Hold the Remote Control as shown above to establish contact with both sound processors

- 1 Searching for two sound processors
- 2 Searching for the right sound processor
Connected to the left sound processor

Display

To conserve battery life, the Remote Control has a power-saving feature that automatically turns off the display and puts the Remote Control in standby mode after 20 seconds of inactivity. Simply press any button on the keypad to reactivate the display and then perform the desired action.

Remote control performance

The Remote Control performs best when buttons are pressed one step at a time rather than in rapid sequence. When adjusting volume, for example, briefly press a button as described below and wait for the tone before continuing. This tip applies to all button use but is especially relevant to volume adjustments.

Press any button on the Remote Control to reactivate it from standby mode. Press again to make the intended changes (e.g., increase volume, change programs, etc.).

Pairing

See figure **10**

Pairing is the process of establishing communication between your sound processor and your Cochlear Baha Remote Control. Your hearing care professional can help you with this procedure when fitting your sound processor. Alternatively, you can follow the steps below to pair your sound processor with the Cochlear Baha Remote Control:

1. Make sure the Cochlear Baha Remote Control is fully charged and is turned on.
2. Press the pairing button on the Cochlear Baha Remote Control once. This will activate the pairing mode **1** which will be active for 20 seconds.
3. While pairing mode is activated, turn on your sound processor. Successful pairing will be indicated by an audible melody played in your sound processor and program 1 will be shown as active in the remote control display **2**. You are now ready to use the Remote Control.

Cochlear TV Streamer (optional)

See figure **11**

The Cochlear TV Streamer is an optional audio and speech streaming device available for your sound processor.

The Cochlear TV Streamer streams audio signals (stereo) from your TV, stereo, computer, etc., directly to your sound processor.

Activating streaming from the Cochlear TV Streamer

See figure **12**

To start audio streaming from the TV, stereo or any other device connected to your Cochlear TV Streamer, simply turn on the audio device and then press the Cochlear TV Streamer button (no. 4 in figure 1). The Remote Control will immediately start searching for the Cochlear TV Streamer and connect to it.

- 1 Searching for a TV
- 2 Connected to a TV
- 3 No TV found

You can have up to 3 Cochlear streaming devices set up for your sound processor. Each Cochlear streaming device is indicated in the display by the letter A, B or C. If more than one Cochlear streaming device is in use, press the streaming button again until the desired audio signal is heard.

To stop streaming and return to the previous program, press the "P" button (no. 6 in figure 1). Alternatively, press the Home button (no. 4 in figure 1), which will also reset the Remote Control to simultaneous adjustment mode. Please note, however, that the Home button also changes other settings, as described above.

Appendix: exceptions to normal operation

If you have been fitted with two sound processors, it is recommended that both of them be set to the same program at all times. However, it is possible to have different programs on the two sound processors, such as Program 1 on the left side and streamed audio on the right. Under these circumstances, the remote control functions slightly differently, as the following examples illustrate.

Situation #1: You have started with P1 on both sides but have changed the right side to P2 by pushing the button on the right sound processor. If you then press the P button on the remote control, the left sound processor will advance to P2, matching the one on the right. Further presses of the P button will advance both sound processors, thus keeping the same program on both sides.

Situation #2: You have set the left sound processor to P1 and the right sound processor to P2. If you then press the Streaming button, both sound processors will receive streamed audio from the same source (e.g., the TV). If you then press the P button on the remote control, streaming will be cancelled on both sides and both sound processors will return to their previous programs (P1 and P2).

Situation #3: You have P1 on the left side but have started streaming on the right side by pressing and holding the button on the right sound processor. If you then press the P button on the Remote Control, streaming on the right side will be cancelled and both sound processors will return to P1.

If you press the Streaming button instead, both sound processors will be set to the same streamed audio program.

Important information

Maintenance and Care

- Clean the Cochlear Baha Remote Control using a damp cloth. Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the device.
- When the Cochlear Baha Remote Control is not in use, turn it off and store it safely.
- Protect the Cochlear Baha Remote Control from excessive moisture

(bathing, swimming) and heat (radiator, car dashboard).

- Protect it from excessive shock and vibration.

General warnings

- Keep this device out of reach of children.
- The Cochlear Baha Remote Control uses low-powered digitally coded transmissions in order to communicate with other wireless devices. Although unlikely, nearby electronic devices may be affected. In such cases, move the Cochlear Baha Remote Control away from the affected electronic device.
- When using the Cochlear Baha Remote Control and it is affected by electromagnetic interference, move away from the source of the interference.
- Avoid exposing your product to rain, moisture or other liquids to protect against damage to the product or personal injury.
- Keep the product and any cords and cables away from operating machinery.
- If the product overheats, has been dropped or damaged, has a damaged cord or plug, or has been dropped in a liquid, discontinue use and contact your hearing care professional.
- Dispose of the product according to local standards and regulations.
- Do not use in areas where there is a risk of explosion.

General precautions

- Your sound processor and the remote control were given a unique communication network code during fitting. This ensures that the remote control will not affect sound processors worn by others.
- High-powered electronic equipment, larger electronic installations and metallic structures may significantly reduce the operating range.
- Use only original Cochlear accessories.
- Do not alter or modify the remote control.
- Opening the remote control might damage it.
- The Cochlear Baha Remote Control may only be repaired by an authorized service centre.
- Do not use the Cochlear Baha Remote Control in areas where RF transmission is prohibited, e.g., in planes, etc.

Battery warning information

- Batteries contain dangerous substances and should be disposed of carefully.
- For your safety and to protect the environment, please note:
 1. DO NOT attempt to recharge batteries that are not rechargeable. They may leak or explode.
 2. DO NOT attempt to dispose of batteries by burning them.
 3. DO NOT place batteries in your mouth. Consult a physician immediately if a battery has been swallowed. This can be harmful to your health.
 4. Keep batteries away from pets or anyone who may not be capable of handling them properly.

Technical specifications

Operating range	Up to 3 meters
Batteries	2 x 1,5V AAAA
Wireless	2,4 GHz
Operating temperature	0 to 55° C
Storage temperature	-20 to 60° C

Temperature test, transport and storage information

The product is subject to various tests in temperature and damp heating cycling between -25° C and +70° C according to internal and industry standards.

Wireless Accessory type designations for models included in this User Manual are:

FCC ID QZ3RC-1, IC 8039C-RC1, IC model: RC-1

Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- Changes or modifications can void the user's authority to operate the equipment.

List of countries:

The products are in compliance with the following regulatory requirements:

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Equipment includes RF transmitter.

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