Cochlear Baha Remote Control 2™

Welcome

Congratulations on your purchase of a Cochlear™ Wireless TV Streamer.

Intended use

The Cochlear Wireless TV Streamer is intended to be a wireless connection to a television or other audio device for a compatible Cochlear sound processor.

Indications:

This accessory has no therapeutic indications in itself. It will support the therapeutic indications of the Cochlear wireless sound processor used with this accessory.

Clinical Benefit:

The Cochlear Wireless TV Streamer is intended to further improve audibility when watching television or listening to sound from other audio device.

Key to device

- Micro-USB power socket
- Pairing button
- Mono/stereo switch
- Analogue mini-jack audio input
- Optical (Toslink) audio input
- Coaxial audio input
- Increase volume
- Decrease volume
- Power light (green)
- Activity light (yellow)
- Wall plug
- Micro-USB plug

Introduction

The Cochlear Wireless TV Streamer will provide you with high quality streamed audio from your TV, stereo, personal computer or similar device directly to your sound processors.

Please read this manual carefully to fully benefit from the Cochlear Wireless TV Streamer. Ask your hearing care professional if you have any questions. For further information and a full list of compatible sound processors you can also visit the www.cochlear.com/compatibility web site.

For your safety, carefully read the chapter "Important information".

Warranty

The warranty does not cover defects or damage arising from, associated with, or related to the use of this product with any non-Cochlear processing unit and/or any non-Cochlear implant. See "Cochlear Global Limited Warranty card" for more details.

Key to symbols

The following symbols will be used throughout this document.

	Caution
	CE-mark
	UK Conformity mark
	Serial Number

1. Make sure that the TV Streamer is turned on and that your sound processors are receiving the streamed audio.

2. Press and hold the pairing button [no. 2 in figure 1] for 3 seconds. The TV Streamer will then enter the adjustable delay mode where the green indicator light on the front of the TV Streamer will flash once every 2 seconds. The yellow indicator light will be on whenever an audio signal is present.

3. Use the "u" and "-" key [no. 7 and 8 in figure 1] on the top of the TV Streamer to increase or decrease the delay of the streamed signal in 10-m steps. The streamed signal can be delayed in the range of 0-250 ms.

4. When the optimal alignment between the audio streamed to the sound processors and the sound coming from the TV speakers or alignment between the streamed audio and the TV picture has been reached, press and hold the pairing button for 3 seconds. Now the TV Streamer will resume normal operation, indicated by the green indicator light being constantly on.

5. The offset set in adjustable delay mode will be stored in the TV Streamer. If you do not exit the adjustable delay mode by pressing and holding the pairing button for 3 seconds, the TV Streamer will automatically resume normal operation after 30 minutes. To reset the delay to 0 ms (factory setting), press and hold the pairing button while simultaneously powering on the TV Streamer.

Pairing the TV Streamer in channel 2

1. Press the pairing button twice. The yellow indicator light on the front of the TV Streamer will now flash twice every 2 seconds to indicate that it is ready to be paired in channel 2.

2. While pairing mode is active, turn on your sound processor. An audible melody or a flashing light (depending on the type of sound processor) will be followed by streamed audio via the TV Streamer.

Pairing the TV Streamer in channel 3

1. Press the pairing button three times. The yellow indicator light on the front of the TV Streamer will now flash three times every 2 seconds to indicate that it is ready to be paired in channel 3.

2. While pairing mode is active, turn on your sound processor.

Streaming

Turn on the TV, stereo, computer or other audio device connected to the TV Streamer. Make sure your sound processors are turned on as well. Streaming may be started in one of two ways:

Option A
1. Press and hold the wireless streaming button on your sound processor for about three seconds.

2. An audible melody or a flashing light (depending on the type of sound processor) will be followed by streamed audio via the TV Streamer.

Option B
If you have a Cochlear Baha Remote Control** (optional), simply press the streaming button [no. 3 in figure 7] on this Baha Remote Control**. For more information, refer to the Baha Remote Control** manual.

Start streaming on a 2nd or 3rd TV Streamer unit

If you have a Cochlear Baha Remote Control** (optional), simply press the streaming button [no. 3 in figure 7] on this Baha Remote Control**. For more information, refer to the Baha Remote Control** manual.

Alternative analogue audio cable installations:

If it is not possible to connect the TV Streamer to the AUDIO OUT/OUTPUT of the audio source using the red and white plugs, the following installation alternatives can be used:

Mini-jack adapter:

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer.

2. Insert the male stereo mini-jack on the audio cable into the female stereo mini-jack input on the TV Streamer [no. 4 in figure 1].

3. Connect the stereo mini-jack adapter to the red and white plugs at the other end of the audio cable.

4. Plug the stereo mini-jack adapter into the audio source's "HEADPHONE" output terminal or alternatively into your computer's audio terminal.

Note: The HEADPHONE output socket is usually located on the back of the audio source. When using the HEADPHONE socket, usually only the person wearing the sound processors will hear the sound and audio source speakers are muted. If muted speakers are not desired, please refer to the first installation, or the one in the following chapter. Your audio source may support speakers and HEADPHONE being on simultaneously; check your audio source manual for instructions.

Scart adapter:

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer.

2. Insert the male stereo mini-jack on the audio cable into the female stereo mini-jack input on the TV Streamer [no. 4 in figure 1].

3. Connect the scart adapter to the red and white plugs at the other end of the audio cable.

4. Plug the scart adapter into the audio source's scart terminal. **Note:** The scart output terminal is usually located on the back or side of the audio source.

Please consult your audio source's manual if additional scart output configuration is needed.

Pairing with your sound processors

Select a separate channel for each device paired with the sound processor.

Pairing flashing pattern	LED blinking pattern	Channel
1		1
2		2
3		3

Note: The TV Streamer supports both analogue and digital audio input. Please read the following sections for information about supported audio formats and how to connect the TV Streamer to your sound processors.

Connecting the TV Streamer to a digital audio source (e.g. TV or stereo)

Note: The TV Streamer supports both analogue and digital audio input. When using digital audio input the following formats are supported:

1. Make sure that the TV Streamer is turned on and that your sound processors are receiving the streamed audio.

2. Press and hold the pairing button [no. 2 in figure 1] for 3 seconds. The TV Streamer will then enter the adjustable delay mode where the green indicator light on the front of the TV Streamer will flash once every 2 seconds. The yellow indicator light will be on whenever an audio signal is present.

3. Use the "u" and "-" key [no. 7 and 8 in figure 1] on the top of the TV Streamer to increase or decrease the delay of the streamed signal in 10-m steps. The streamed signal can be delayed in the range of 0-250 ms.

4. When the optimal alignment between the audio streamed to the sound processors and the sound coming from the TV speakers or alignment between the streamed audio and the TV picture has been reached, press and hold the pairing button for 3 seconds. Now the TV Streamer will resume normal operation, indicated by the green indicator light being constantly on.

5. The offset set in adjustable delay mode will be stored in the TV Streamer. If you do not exit the adjustable delay mode by pressing and holding the pairing button for 3 seconds, the TV Streamer will automatically resume normal operation after 30 minutes. To reset the delay to 0 ms (factory setting), press and hold the pairing button while simultaneously powering on the TV Streamer.

Important information

Maintenance and Care

• Clean the TV Streamer using a damp cloth. Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the device.

• When the TV Streamer is not being used for extended periods, unplug the power supply from the wall socket.

General warnings

• Keep this device out of reach of children.

• The TV Streamer uses low-powered digitally coded transmissions in order to communicate with other wireless devices. Although unlikely, nearby electronic devices may be affected. In such cases, move the TV Streamer away from the affected electronic device.

• When using the TV Streamer and it is affected by electromagnetic interference, move away from the source of the interference.

Power supply warning

• Use the TV Streamer only with the power supply unit that was supplied with the product.

• The use of other power suppliers may damage or destroy the product and could be dangerous. Use of other power suppliers may invalidate any approval or void the warranty.

• When you disconnect the power cord or any accessory, grasp and pull the plug, not the cord.

• Never use a power supply that is damaged.

• Do not attempt to disassemble the power supply because it may expose you to a dangerous electric shock.

• Do not use the power supply outdoors or in damp areas.

General precautions

• Do not use this device in locations where the use electronic devices are prohibited.

• X-ray radiation (e.g. CT/MRI scan) may adversely affect the functioning of this device. Switch it off before undergoing the procedure and keep it outside the room.

• Your sound processors and this device were given unique communication codes during fitting. This ensures that the device will not affect sound processors worn by others.

• High-powered electronic equipment, larger electronic installations and metallic structures may impair and significantly reduce the operating range.

• Use only original Cochlear accessories.

• Do not alter or modify this device.

• Opening the TV Streamer might damage it. If problems occur that cannot be resolved, consult your hearing care professional.

• The TV Streamer may only be repaired by an authorized service centre.

• Do not use the TV Streamer in areas where RF transmission is prohibited, e.g. in planes, etc.

• Only connect the TV Streamer to connections for which it is explicitly intended.

Serious incidents

Serious incidents are rare, any serious incident in relation to your device should be reported to your Cochlear representative and to the medical device authority in your country, if available.

Technical specifications

Operating range	Up to 7 meters
Power supply	External wall adapter
Wireless	2,4 GHz
Operating temperature	0 to 55 °C
Storage temperature	-20 to 60 °C

Temperature test, transport and storage information

The product is subject to various tests in temperature and damp heating cycling between -25 °C and +70 °C according to internal and industry standards.

Wireless Accessory type designations for models included in this User Manual are:

Model: SAS-3
 FCC ID: QZ3SAS-3
 IC: 8039C-SAS3

 CCAH15LP1740T5
 CMIIT ID:2015DJ4572

Adjustable delay mode

In rare cases there may be a misalignment between the audio streamed to the sound processors and the sound coming from the TV speakers (echo) or even a misalignment between the streamed audio and the TV picture (lip-synch). This may be evident in more complex audio setups where the TV Streamer is not connected directly to the TV (e.g. through a set-top box or other auxiliary unit). To overcome such misalignments the streaming delay on the TV Streamer can be adjusted as follows:

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