

Patient Identification Card

This person is implanted with a medical hearing device manufactured by Cochlear Limited. This device is made of metal and implanted under the skin behind the ear. It usually contains a magnet and is connected to electrodes. External components are worn behind the ear, on the head or on the body, e.g. in a pocket. Without these components, this person cannot hear.

These components may activate metal detection devices, e.g. airport security detector.



The patient must remove all external components of their Cochlear implant system before entering a room where an MRI scanner is located.

The implant is MR Conditional. The implant model information is available on the reverse side of this card. MRI examinations can be performed safely on this person only under very specific conditions. MRI examinations performed under different conditions may result in severe patient injury or device malfunction.

Full MRI safety information for this implant is available by:

- visiting www.cochlear.com/warnings
- calling your regional Cochlear office – contact numbers are available on the reverse side of this card

Hear now. And always



First name: _____ Last name: _____

Address: _____

_____ Country: _____ Tel: _____

Emergency contact: _____ Tel: _____

Implant serial no.: _____ Model: _____ Ear: _____

Implant serial no.: _____ Model: _____ Ear: _____



Cochlear Ltd (ABN 96 002 618 073), Australia, Tel: +61 2 9428 6555

Cochlear Americas, USA, Recipient Customer Service, Tel: (800) 483 3123

Cochlear Canada Inc, Canada, Tel: +1 416 972 5082

www.cochlear.com

537275 ISS3