




















Cochlear™ Baha® compatibility matrix: Sound processors, Wireless accessories, Apps and MFi












	Mini Microphones			TV Streamer	Remote Controls		Phone Clip	Apps		MFi
										
	Mini Microphone 1	Mini Microphone 2	Mini Microphone 2+	TV Streamer	Remote Control 1	Remote Control 2	Phone Clip	Control App	Smart App	Made for iPhone
 Baha 5 Sound Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
 Baha 5 Power	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
 Baha 5 SuperPower	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
 Baha 4 Sound Processor	✓	✓	✓	✓	✓	✓	✓	✓		
Operating system requirements							iOS 7.X or later.	iOS 7.X or later. Android OS 2.3 or later.	iOS 9.1 or later Android OS 5.0 (Lollipop) or later	iOS 7.X or later.
Hardware requirements							iPhone X, iPhone 8 Plus, iPhone 8, iPhone 7 Plus, iPhone 7, iPhone 6s Plus, iPhone 6s, iPhone 6, iPhone 6 Plus, iPhone 5s, iPhone 5c, iPhone 5, iPhone 4s and iPhone 4.	iPhone, iPad or iPod touch and Android smartphones.	iPhone, iPad and iPod touch Samsung Galaxy S6, S7 and S8.	iPhone X, iPhone 8 Plus, iPhone 8, iPhone 7 Plus, iPhone 7, iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 3, iPad mini 2, iPad mini and iPod touch (5th generation).
Other Requirements								Cochlear Wireless Phone Clip.		

Cochlear™ Baha® compatibility matrix: Systems

Baha Connect System

Implants	Abutments	Sound processors
 <p>BI300 3 mm, 4 mm</p>	 <p>BA300 6 mm, 9 mm</p>  <p>BA400 6 mm, 8 mm, 10 mm, 12 mm, 14mm</p>	 Baha 5 Sound Processor  Baha 5 Power  Baha 5 SuperPower  Baha 4 Sound Processor  Baha BP100  Baha BP110 Power  Baha Cordelle II
 <p>Generation 2 3 mm, 4 mm</p>	 <p>Generation 2 5.5 mm, 8.5 mm</p>  <p>BA210 5.5 mm, 8.5 mm</p>	 Baha 5 Sound Processor  Baha 5 Power  Baha 5 SuperPower  Baha 4 Sound Processor  Baha BP100  Baha BP110 Power  Baha Cordelle II

Baha Attract System

Implants	SP Magnets	Sound processors
 <p>BIM400</p> <p>+</p>  <p>BI300 3 mm, 4 mm</p>	 <p>SP Magnet Generation 2 Strengths 1 – 6</p>  <p>SP Magnet Generation 1 Strengths 1 – 6</p>	 Baha 5 Sound Processor  Baha 5 Power  Baha 5 SuperPower  Baha 4 Sound Processor  Baha BP100  Baha BP110 Power  Baha Cordelle II

www.cochlear.com

Baha, Cochlear, Hear now. And always and the elliptical logo are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB or Cochlear Limited. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. Android is a trademark of Google Inc. iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. Samsung, Galaxy S6, Galaxy S7 and Galaxy S8 are trademarks of Samsung Electronics Co. Ltd. © Cochlear Bone Anchored Solutions AB 2018. All rights reserved. FEB18.